

VLOSTIC LIGHT hair tonic

VLOSTIC LIGHT Is a light tonic for hair and scalp care. Thanks to the multifunctional formula it prevents hair loss and strengthens their structure. The preparation VLOSTIC LIGHT detangles the hair, gives them shine and smoothness and makes it easier to comb.

INSTRUCTIONS OF USE

DAILY CARE FOR HAIR AND SCALP. VLOSTIC LIGHT IS INTENDED FOR EVERYDAY CARE OF:

- ✓ **DAMAGED HAIR,**
- ✓ **HAIR WEAKENED BY COSMETIC AND MEDICAL PROCEDURES,**
- ✓ **HAIR WITH TENDENCY TO LOSS.**

AN UNIQUE FORMULA

- ✓ **VITAMIN COMPLEX** containing vitamins A, E, H, B5, B8 and vegetable oils rich in omega-6 and omega-3.
- ✓ **EXTRACT OF ROSEMARY** strengthens hair roots, protects the hair fiber from damage and fading, and gives them shine and smoothness.
- ✓ **EXTRACT OF ARNICA MONTANA** stimulates blood circulation in the scalp, protects the blood vessels.
- ✓ **EXTRACT OF BURDOCK ROOT** has a strong bacteriostatic effect.
- ✓ **EXTRACT OF GINSENG ROOT** stimulates hair growth.
- ✓ **EXTRACT OF EQUISETUM ARVENSE** strengthens and revitalizes hair structure.
- ✓ **EXTRACT OF IVY LEAFS AND SHOOTS** improves microcirculation and strengthens blood vessels.
- ✓ **EXTRACT OF WATERCRESS LEAFS AND FLOWERS** has a bacteriostatic effect.
- ✓ **PANTHENOL** moisturizes, smoothes and thickens the hair fiber
- ✓ **CONDITIONING COMPLEX** provides optimal smoothness and softness, gives an immediate effect without weighting down.

HOW TO USE

Apply the product to wet or to dry hair, gently massage and comb. For daily use.

SIDE EFFECTS

Pruritus, erythema, swelling of the skin or subcutaneous tissue. In the case of occurrence of any of the symptoms one should immediately stop using the product and contact a doctor.

CONTRAINDICATIONS

Hypersensitivity to tonic ingredients, pregnancy, lactation, age under 3 years. Do not use on damaged skin, mucous membranes, genitourinary area, eyes and skin around the eyes.

RECOMMENDATIONS AND PRECAUTIONS

For external use only. If swallowed, immediately contact a doctor. Keep out of the reach of children. Do not use after the expiration date, printed on the package.

Store in a dry place, at a room temperature (15 ° C, - 25 ° C.).

Once opened, use within 6 months.

Dermatologically tested.

Available packaging: bottle with a capacity of 100 ml.

LEFROSCH SP. Z O.O.

ŚNIATY, WIELICHOWSKA STREET 29 64-050 WIELICHOWO

SALES DEPARTMENT JEZIORNA STREET 1 62-052 KOMORNIKI

tel. +48 61 893 53 83,

e-mail: lefrosch@lefrosch.com

www.lefrosch.com